

Atomic Absorption Spectrometer AAS6000



Atomic Absorption Spectrometer

AAS6000

Product Introduction



A New Standard of Healthy Life



AAS6000 Series are Single Beam Atomic Absorption Spectrometers controlled and data processed by external computer and internal CPU chips. AAS6000 can be used to measure the concentration of macro, micro and trace metal elements and half-metal elements in various kinds of substances.

Atomic absorption spectrometry is a technique exploiting the fact that atoms at ground state may absorb the radiation of characteristic wavelength. As is known, the atoms are usually at ground state. For each element, the amount of energy required for an atom to transit from ground state to excited state is specific, and in general, is called characteristic wavelength. The Atomic absorption spectrometry typically makes use of hollow cathode lamp as the light source to emit the characteristic wavelength of a particular element. When the light passes through the atomic vapor, it is absorbed by the vapor. The concentration of the element can then be determined through the calculation of the absorption rate.

AAS6000 is provided with three reading methods Continuum, Retention and Peak Height for measuring absorbency, density and emission intensity. It has three signal modes: Atomic Absorption, Background Absorption and Background Correction Absorption. The reading time ranges from 0.5 s to 99s. AAS6000 include 8 computer controllable hollow cathode lamp holders. All working conditions of the instrument such as lamp number, lamp current, negative high voltage, working wavelength, slit, burner location, ignition/extinction and adjustment of burning ratio can be set by keyboard input. The functions of the instrument include automatic gain/ automatic zero, background correction, automatic energy balance, automatic peak location/wavelength scan, and automatic peak location on basis of the retrieval of peak values. All the readings, measurement results, calibration curves and operation conditions can be saved or printed out.

Application fields

- Metallurgy, steel, non-ferrous metals.
- Environmental analyses: air, water quality, soil and solid wastes.
- Petrochemical engineering, crude oil and related products, light industrial products.
- Food, biomedicine and health products.
- Building materials (glass, ceramic, paints, etc).



Atomic Absorption Spectrometer

AAS6000

Atomic Absorption Spectrometer AAS6000



Periphery equipments

PC (Brand computer, P4, LCD)	1 set
Printer (Canon, colorful ink-jet printer)	1 set
Air compressor	1 set

Main parts

Narrow-line light source
Gas control system
Optical system
Circuit control system

Optional components

AC stabilized power supply
Ultrasonic nebulizer
Hydride generator

Power requirements

Single phase 220V±10%, 10A
One switchboard (220V, 10A) is needed. The instrument should be within the 10m range of the switchboard, and earthed.
Two plug boards (220V, 10A)

Gas

Purity of acetylene
Navigation mark acetylene
Usage of acetylene
1. 6L/min@0.2MPa

Specifications

Model: AAS6000
Optical system: all-reflective single-beam CT optical path
Focal length of monochromator: 350mm, oil/water proof air compressor
Blazed wavelength of gratings: 230nm
Number of grating grooves: 1800g/mm
Wavelength range: 190nm-900nm
Wavelength accuracy: ±0.1nm
Wavelength repeatability: ±0.1nm
Spectrum bandwidth: 0.1/0.2nm/0.4nm/0.7nm/1.4nm
Noise: 0.005 Abs(Static); 0.006 Abs(Dynamic)
Baseline drift: 0.003 Abs/0.5h; best performance
Background correction: D₂ lamp + self reversal
Number of the lamps: 8
Number of preload lamps: software controlled, ≤8
Gas path safety measures: yes
Atomizer: flame atomizer
Automation: lamp/slit/wavelength/ gas path/ignition/burner/protection/ C₂H₂ monitoring

Advanced features

Advanced features:

- C₂H₂ monitoring; flame measurement; all reflective, aberration corrected.
- Fully automated instrument, total automatic operation.
- 8 element lamps, time saved for changing lamps.
- Gas path safety measures, safer operation.
- USB interface allows you to operate the instrument conveniently and rapidly.
- 230nm grating blazed wavelength, multiplied sensitivity in the Ultra violet zone.
- 1800g/mm grating grooves, maximized resolution.
- Run in Windows operating system, the powerful software enables automatic qualitative and quantitative analysis. Automatic report generation in Excel format allows you to operate the instrument and process the data easily. Your operation does get any easier with our graphical interfaces.



Test example

The following Figure is a test example of Pb content in liquid solvent taken with Skyray self- developed AAS6000 Atomic Absorption Spectrometer. With figurative interfaces and total automatic measurement, the instrument offers best test results. The advanced features of AAS6000 include: high sensitivity, low detection limit, good selectivity, uncompromised accuracy, simple operation, elemental range of 68 elements including both metals and non metals and organic chemicals.








Figure 1 is good representation of the performance of AAS6000: good baseline stability, accurate wavelength peak location, powerful background correction, rapid, high efficient and precise analysis of liquid solvent, superior test results, intuitive display of instrument parameters, real-time absorbency and work curve, providing a better way to view the measurement results.

Figure 2 is energy spectrum, showing advanced features of our product as well: accurate wavelength, good repeatability, high resolution, intuitive display of instrument parameters and so on.

Note: As the hardware and software are in debug phase, the specifications of the instrument may be subject to changes without notice. Your kind understanding is highly appreciated.



Test report

	Element: Pb Method: flame absorption		Wavelength(nm): 283.30																																																																										
Instrument parameters	Spectral bandwidth (nm): 0.2 Lamp current (mA): 3.00 Flame type: Nitrous oxide -acetylene Fuel gas pressure (MPa) :0.10 Support gas flow rate (L/min):5.00 Height of burner head (mm):20.0		Negative high voltage (V) :256 Background correction Fuel gas flow rate (L/min):1.60 Support gas: air Support gas pressure (MPa):0.20																																																																										
Analysis parameters	Sampling speed: 50 Delay(s):0 Ruler extension (0.1-100) :1.00		Integral time(s) :1 Calculation method: peak height Flame micro sample introduction: no																																																																										
Spectrum	<div><div></div><div></div></div> <div>Linear correlated coefficient: 0.9996 Curve equation: Y=0.0166*X+0.0082</div>																																																																												
Analysis results	<table><thead><tr><th>Name</th><th>Times</th><th>ABS</th><th>Density (ppm)</th><th>SD</th><th>RSD (%)</th><th>Results*</th></tr></thead><tbody><tr><td rowspan="4">Standard 1</td><td>1</td><td>0.0074</td><td rowspan="4">0.0000</td><td rowspan="4">0.00052</td><td rowspan="4">7.50411</td><td rowspan="4"></td></tr><tr><td>2</td><td>0.0069</td></tr><tr><td></td><td>0.0064</td></tr><tr><td>3</td><td>0.0069</td></tr><tr><td rowspan="4">Standard 2</td><td>1</td><td>0.0949</td><td rowspan="4">5.0000</td><td rowspan="4">0.00078</td><td rowspan="4">0.082395</td><td rowspan="4"></td></tr><tr><td>2</td><td>0.0938</td></tr><tr><td></td><td>0.0935</td></tr><tr><td>3</td><td>0.0941</td></tr><tr><td rowspan="4">Standard 3</td><td>1</td><td>0.1748</td><td rowspan="4">10.0000</td><td rowspan="4">0.00078</td><td rowspan="4">0.76521</td><td rowspan="4"></td></tr><tr><td>2</td><td>0.1722</td></tr><tr><td></td><td>0.1730</td></tr><tr><td>3</td><td>0.1733</td></tr><tr><td rowspan="4">Sample 1</td><td>1</td><td>0.0962</td><td rowspan="4">5.2470</td><td rowspan="4">0.00087</td><td rowspan="4">0.90879</td><td rowspan="4">5.25</td></tr><tr><td>2</td><td>0.0945</td></tr><tr><td></td><td>0.0950</td></tr><tr><td>3</td><td>0.0953</td></tr><tr><td rowspan="4">Sample 2</td><td>1</td><td>0.0083</td><td rowspan="4">0.0723</td><td rowspan="4">0.00091</td><td rowspan="4">9.66425</td><td rowspan="4">0.07</td></tr><tr><td>2</td><td>0.0097</td></tr><tr><td></td><td>0.0101</td></tr><tr><td>3</td><td>0.0094</td></tr></tbody></table>					Name	Times	ABS	Density (ppm)	SD	RSD (%)	Results*	Standard 1	1	0.0074	0.0000	0.00052	7.50411		2	0.0069		0.0064	3	0.0069	Standard 2	1	0.0949	5.0000	0.00078	0.082395		2	0.0938		0.0935	3	0.0941	Standard 3	1	0.1748	10.0000	0.00078	0.76521		2	0.1722		0.1730	3	0.1733	Sample 1	1	0.0962	5.2470	0.00087	0.90879	5.25	2	0.0945		0.0950	3	0.0953	Sample 2	1	0.0083	0.0723	0.00091	9.66425	0.07	2	0.0097		0.0101	3	0.0094
Name	Times	ABS	Density (ppm)	SD	RSD (%)	Results*																																																																							
Standard 1	1	0.0074	0.0000	0.00052	7.50411																																																																								
	2	0.0069																																																																											
		0.0064																																																																											
	3	0.0069																																																																											
Standard 2	1	0.0949	5.0000	0.00078	0.082395																																																																								
	2	0.0938																																																																											
		0.0935																																																																											
	3	0.0941																																																																											
Standard 3	1	0.1748	10.0000	0.00078	0.76521																																																																								
	2	0.1722																																																																											
		0.1730																																																																											
	3	0.1733																																																																											
Sample 1	1	0.0962	5.2470	0.00087	0.90879	5.25																																																																							
	2	0.0945																																																																											
		0.0950																																																																											
	3	0.0953																																																																											
Sample 2	1	0.0083	0.0723	0.00091	9.66425	0.07																																																																							
	2	0.0097																																																																											
		0.0101																																																																											
	3	0.0094																																																																											